

**Product:** CEDERROTH EYE WASH  
POCKET MODEL  
**Reference code:** 7221



## 1. General description

A personal eyewash bottle containing sterile, buffered, saline eye wash solution.

The bottle is easily opened and ready to use by twisting off the dustcover. The integrated eye cup helps keeping the eye open during rinsing and ensures a generous flow without unnecessary waste. The bottle is shaped to fit into a pocket.

## 2. Intended use

The Cederroth Eye Wash Pocket Model is intended for washing away foreign particles and dangerous liquids, like alkali or acids, from an eye in an emergency first aid situation at work places.

The Cederroth buffered Eye Wash also has a neutralising effect on splashes from alkalis and acids. The effect is in general stronger for alkalis than for acids.

The Cederroth Eye Wash is intended for the first acute phase at the very accident site before getting to a plumbed or self-contained eye wash equipment or seeking medical attention.

## 3. Technical description

### *Product dimensions :*

(l x w x h) mm: 45 x 66 x 163

Fill volume (ml): 235

### *Material*

Liquid Sterile, isotonic Sodium Chloride 0.9%,  
Buffered with Sodium Borate to neutral PH.  
Preservative free.

Bottle Leak proof polypropylene bottle produced with Blow-Fill-Seal technique.  
Eye cup in transparent polypropylene  
Dust cover in green hard plastics (ABS)

## 4. Sterilization

Sterilization method: Aseptic filling

## 5. Composition

### *Consumer package*

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>
7221	Cederroth Eye Wash Pocket Model	1

### *Transport package*

Number of consumer packages in transport package:	18 (3 x 6 ) pcs
Dimensions (l x w x h) mm:	477 x 162 x 173
Weight (kg):	5,4

## 6. Labeling

### *On product/ consumer package*

Cederroth address, telephone- and fax-number, e-mail and web-address  
 Product name, Reference code  
 Lot-number, Expiry date, EAN-code  
 Instructions for use  
 Symbols: CE +NB#, STERILE A,  
 DO NOT REUSE, LOT, USE BY,  
 TEMPERATURE LIMITATION

### *On transport package*

Cederroth address, telephone- and faxnumber,  
 web-address  
 Product name, Reference code  
 Lot-number, Expiry date, EAN-code  
 Quantity  
 Symbols: LOT, USE BY

## 7. Storage conditions

Store in temperatures between +2°C and +40°C.

Shelf life: 4,5 years from production date

## 8. Accessories

The following accessories can be used with the product but are not necessary.

720300

Holster for Cederroth Eye wash pocket model



## Authorization: