

PRODUCT SPECIFICATION



Product: CEDERROTH 4-IN-1

Reference code: 1910



1. General description

Cederroth 4-in-1 is a sterile universal dressing, which consists of two elastic bandages and one compress with a non-stick wound contact layer. The compress is stitched to one of the bandages. The dressing is wrapped in a pouch with printed pictorial user instructions.

2. Intended use

To be used in first aid situations. Cederroth 4-in-1 has four functions:

- Pressure pad for heavy bleeding
- > Protective dressing for superficial wounds
- ➤ Burn dressing with special wound protecting surface
- > Temporary support dressing

Because of the many functions, Cederroth 4-in-1 replaces a variety of traditional first aid dressings and compresses.

3. Technical description

Product dimensions:

Consumer package: 12,5 x 6,5 x 3 cm Bandage: 300 x 10 cm Compress: 23 x 14 cm

Material

Bandages Rayon

Compress Absorbent inner layer: Rayon staple fibers

Non adherent outer cover: Porous net made of high density

polyethylene

Packing material Transparent film made of polyethylene. Tyvek.

(sterile package)

4. Sterilization

Sterilization method: Irradiation

5. Composition

Consumer package Transport package

Ref#	Product/Material	<u>Pcs</u>	Number of consumer packages
51616500	Big bloodstopper	1	in transport package:

51616515 Plastic pouch 1 Dimensions (1 x w x h) cm: 26 x 21 x 18

Weight (kg):

20



PRODUCT SPECIFICATION



6. Labeling

On product/ consumer package

Cederroth address, telephone- and faxnumber, web-address, mail-address

Product name, Reference code Lot-number, Expiry date, EAN-code, Single use, Sterile R, User instructions, Product size, CE+NB#

On transport package

Cederroth address, telephone- and faxnumber, web-address
Product name, Reference code
Lot-number, Expiry date, EAN-code, Quantity, CE+NB#

7. Storage conditions

The product should be stored in the original package for maximum protection until time of use. Do not use if the sterile package is broken.

Shelf life: 5 years from production date

Authorization: